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FAUTAPO

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Cocoa Research Centre
The University of the West Indies
St. Augustine, 330912, Trinidad and Tobago
www.cocoacentre.com; www.ifcic.com

Tel/fax +1 868 662 8788

E-mail: lambert.motilal@sta.uwi.edu

Sending Samples in Collection Tubes.


1. Collection tubes contain desiccant beads which are yellow or golden yellow in colour. Please examine collection tubes on arrival. If the beads appear light yellow – colourless, remove beads from tubes and heat in microwave or oven (not more than 100°C) until the golden yellow colour is restored. Let beads cool before returning to tubes.
2. Keep tubes sealed and away from moisture until use.
3. Each collection tube should contain samples from one tree, for example leaf pieces from one leaf of one plant. Write a unique code on each collection tube to match that of each unique sampled plant.
4. Examine the leaf for a particular tube and avoid any yellow, brown or diseased areas.
5. Have a clean pair of scissors, alcohol and paper towels available. Clean the scissors blade with an alcohol-soaked paper.
6. Wipe both upper and lower leaf surfaces with alcohol.
7. Cut panels of leaf tissue between the major veins of the leaf as in the figure below:

8. Size the leaf panels to fit the collection tube.

9. Gently curve the sections so that they can fit easily as a single layer and without much overlap within the collection tube. Place curved leaf section on top of the desiccant (orange beads). If the collection tube is large enough, you may put more than one panel inside the collection tube. However, a collection tube must contain samples from only one tree.

10. Close collection tube securely. The panel(s) may rest on the desiccant beads or be shaken to distribute the beads around the panel(s).
11. Using a permanent marker, write a unique label on the collection tube to match that of the sampled tree.
12. Clean the scissors with alcohol before moving on to the next sample. Ensure that the blades are free of any liquid or debris from the cut leaf of the previous sample.
13. Keep labelled sample collection tubes away from direct sunlight and moisture.
14. Double bag the sample collection tubes in zipper bags. Remove excess air from bags.
15. Place bagged sample tubes and a packing list (number, names & location of samples) in a box with packaging material.
16. Place separate labels on outside of box with (a) sender, (b) recipient contact details and (c) statement that package is for research purposes only.
17. Attach import permit and phytosanitary certificate as needed and ship to the address above.

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